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

FELDMAN AORTIC STENOSIS CATHETERS

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgement concerning patient care.

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Feldman Aortic Stenosis Catheter is available in three curve sizes, -FA, -FA1 and -FA2 (Fig. 1). The nylon catheter with stainless steel braiding (HNB-) provides excellent memory and torque control to the catheter tip. The catheter tip has two opposing sideports and a tapered end hole for use during ventriculography.

The FA is used for small aortic roots; FA1 for normal sized aortic roots; FA2 for extremely large aortic roots or extremely angulated valves. The FA1 curve configuration is appropriate for most patients.

Refer to label for further information.

2 ENGLISH

I-FA-0312-224-02

INTENDED USE

The product has been designed for ventriculography. The catheter shape design, coupled with the use of a straight movable core wire guide, facilitates rapid crossing of stenotic aortic valves.

The products are intended for use by physicians trained and experienced in catheterization techniques and angiographic procedures.

Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

CONTRAINDICATIONS

None known.

WARNINGS AND PRECAUTIONS

- Possible allergic reactions should be considered.
- Manipulation of products requires fluoroscopic control.
- Do not exceed recommended flow rate. Flow rates in excess of 10 cc/sec have been reported to cause ventricle perforation and intramuscular injection with improperly positioned catheters. Recommended flow rates of 6 cc/sec for a total volume of 25 cc to 35 cc is sufficient for most patients.
- Carefully position catheter tip in the mid cavity of the left ventricle. Obtain images in the RAO and LAO views to ensure proper positioning toward the apex of the ventricle. A test injection of contrast at a rate of 6 cc/sec for a total volume of 5 cc to 10 cc is recommended to prevent staining of the left ventricle prior to power injection.
- Insertion through synthetic grafts or severely scarred areas should be avoided. If necessary, use of an introducer sheath is recommended.
- Do not attempt to heat or reshape the catheter curve; the catheter tip is made from a heat sensitive material.

RECOMMENDED FLOW RATES

Use of 6 cc/sec for a total volume of 25 cc to 35 cc is adequate for most patients. When the ventricle is extremely dilated, 10 cc/sec may be needed.

I-FA-0312-224-02

ENGLISH 3

INSTRUCTIONS FOR USE

1. Insert the catheter into the aortic root with a straight wire guide and rotate clockwise to direct the tip of the catheter toward the center of the aortic root. No more than 90-120 degrees of rotation is necessary.


2. A nontapered straight movable core wire guide may be used to change the angle of the catheter allowing the operator to scan or interrogate the surface of the aortic valve (Fig. 2). By withdrawing the core of the movable core wire guide 7.0 - 10.0 cm, the tip of the wire guide will be made more flexible, thereby allowing the catheter curve to assume its formed curve completely.

When the aortic root is small in diameter, the catheter can be straightened with the wire and formed into a less-curved shape which approximates a right Judkins curve. The angled design will reach the center of the aortic root in most patients unless the left ventricle and aortic root meet at an extremely acute angle.

4 ENGLISH

I-FA-0312-224-02

3. Once the tip of the catheter and wire guide have passed through the stenotic valve, advance the catheter into the mid-portion of the left ventricle with the tip of the catheter positioned toward the apex. Remove the wire guide. Perform a test injection of 6 cc/sec for a total volume of 5 cc to 10 cc to confirm proper catheter positioning within the left ventricle. Both RAO and LAO projections should be checked prior to power injection. If the catheter tip is positioned against the ventricle wall, reposition catheter tip pointing toward the apex and repeat test injection.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

I-FA-0312-224-02

ENGLISH 5